

Materials Needed to Repair a LDA Using the Grymer / Sterner Toggle Suture Method

- Restraint and Casting Ropes
- * Tranquilizer/Sedation (xylazine/Rompun, BAYER) (xylazine/AnaSed, Lloyd Laboratories)
- Surgical Site Preparation (Scrub and Disinfectant)
- Trocar, Cannula, and Push rod (Jorgensen Laboratories, Inc.) <http://www.jorvet.com>
- 1 pair (2) Toggle Sutures (Jorgensen Laboratories, Inc.) <http://www.jorvet.com>
- 2 needle holders
- Stethoscope

* In the United States, the use of xylazine constitutes extra-label usage, and as such, must be used within the context of a valid Veterinarian Client-Patient Relationship (VCPR) in order to comply with the provisions of the Animal Medicinal Drug Use Clarification Act (AMDUCA).